



General Medical Equipment

ENT Equipment
Suction Devices
Ultrasonic Cleaners

General-Cat/GB-201309 V.1

Sturdy Industrial - Making a Safer, Hygienic Planet

Sturdy reserves the right to improve the products herein without prior notice

www.sturdy.com.tw

Ultrasonic, Ultra-Clean

Sturdy Ultrasonic Cleaners can burst waves up to 46,000 cycles per second, producing longitudinal wave forms that cause a build up of dielectric pressure, which results in micro-vacuum bubbles that explode with forceful cleansing action, easily removing dirt and grime from the most hard-to-reach corners and crevices.



Model	Internal Tank Size (LxWxH)	External Dimensions	Capacity	Tank	Cover	Basket	Housing Material	Ultrasonic Output	Heater Ouput	Electronic Timer
UC-400	495x295x150 mm	530x325x400 mm	22L	St. St. #304	St. St. #304	St.St. #304	Aluminum alloy (powder coating)	400W / 40 KHz	400W	1-99 mins.
UC-300	295x240x150 mm	320x265x280 mm	10.5L	St. St. #304	St. St. #304	St.St. #304	Aluminum alloy (powder coating)	300W / 46 KHz	400W	1-30 mins.
UC-150	295x150x100 mm	320x175x240 mm	4.5L	St. St. #304	St. St. #304	St.St. #304	Aluminum alloy (powder coating)	150W / 46 KHz	200W	1-30 mins.
UC-80	150x135x100 mm	175x160x210 mm	2L	St. St. #304	St. St. #304	St.St. #304	Aluminum alloy (powder coating)	80W / 46 KHz	N/A	1-30 mins.
UC-50	168x89x55 mm	193x116x118 mm	0.6L	St. St. #304	PS Plastic	ABS Plastic ²	ABS Plastic ²	50W / 46 KHz	N/A	N/A

² U/L Fire Resistant, Power Switches on all machines are soft touch buttons, Power Source 110V or 220V or 240V (as requested on order)