

Falcon Series

Vertical Chamber for large instruments

- Pressure Overload Protection
- Overheat Protection
- High mobility



SA-300VF

STURDY
APEX GROUP

Sturdy Industrial, a member of Apex Group, has been supplying reliable sterilization systems for more than 25 years. Sturdy's autoclave is 100% MIT (Made in Taiwan), ensuring every item is designed, developed and manufactured under Total Quality Control standards that surpass strict ISO and CE regulations.



Chamber Capacity	50L
Sterilization Temp. °C	118-134
Dry Time	0-60 Minutes

Other autoclave

- CSSD Solutions
- Lab Solution



With over 360 partners in 60 countries.

Sturdy put safety and quality as first priority.
Trusted in both medical and non-medical industry.

Apex Medical Thailand
65/100 Chamnanphenjati Business Center Bldg.
11/F Rama 9 Road, Huaykwang
Bangkok 10320

Contacto:

Print-2018/07-1000

All rights reserved Specifications may change without notice.

www.amt.in.th

STURDY
APEX GROUP

Reliable Sterilization Partner



Poleax Series

All about the basics

- Overheat Protection
- Pressure Overload Protection
- Pressure Safety Valve

SA-232



Chamber Capacity	16L
Sterilization Temp. °C	126

- Door Safety Latch
- Pressure auto door lock

SA-232X



Chamber Capacity	16L
Sterilization Temp. °C	118-134

Falcon Series

All-time Bestseller

- With Dry Function
- Pressure overload protection
- Overheat protection
- Pressure auto door lock
- Emergency stop button
- Chamber Volume of 40L

SA-300H



Chamber Capacity	40L
Sterilization Temp. °C	118-134
Dry Time	0-60 Minutes

Falcon Series

- Smaller version of SA-300H
- Ideal for little to medium usage

SA-252F



Chamber Capacity	24L
Sterilization Temp. °C	118-134
Dry Time	0-60 Minutes

Super Microm Series

- Pre & Post vacuum
- B.D. & Pump Test Diagnostics
- Prion preset program
- Auto water fill
- With Dry Function
- Optional: Recorder/Flash Cycle

SA-260MA



Chamber Capacity	24L
Sterilization Temp. °C	121 & 134
Dry Time	0/15/30/35/40 Minutes

